



Pragostem

ZYTOVISION
Molecular diagnostics simplified


FISH sondy ZytoLight®

Prenatální, postnatální a preimplantační genetika



FISH sondy ZytoLight®

Prod. No.	Product	Alpha/Class. Sat.	Chr. Band	Label	Tests* (Volume)
Z-2101-200	ZytoLight SPEC 1p12 Probe C€ IVD	-	1p12	●	20 (200 µl)
Z-2049-200	ZytoLight SPEC 2q11 Probe C€ IVD	-	2q11.2	●	20 (200 µl)
Z-2001-200	ZytoLight CEN 3 Probe C€ IVD	D3Z1	3p11-q11.1	●	20 (200 µl)
Z-2083-200	ZytoLight SPEC 4p11 Probe C€ IVD	-	4p11	●	20 (200 µl)
Z-2002-200	ZytoLight CEN 6 Probe C€ IVD	D6Z1	6p11.1-q11.1	●	20 (200 µl)
Z-2003-200	ZytoLight CEN 7 Probe C€ IVD	D7Z1	7q11.1	●	20 (200 µl)
Z-2004-50/-200	ZytoLight CEN 8 Probe C€ IVD	D8Z2	8p11.1-q11.1	●	5/20 (50/200 µl)
Z-2067-200	ZytoLight CEN 9 Probe C€ IVD	III D9Z3	9q12	●	20 (200 µl)
Z-2079-200	ZytoLight CEN 10 Probe C€ IVD	D10Z1	10p11.1-q11.1	●	20 (200 µl)
Z-2005-200	ZytoLight CEN 11 Probe C€ IVD	D11Z1	11p11.1-q11	●	20 (200 µl)
Z-2050-200	ZytoLight CEN 12 Probe C€ IVD	D12Z3	12p11.1-q11	●	20 (200 µl)
Z-2085-200	ZytoLight SPEC 13q12 Probe C€ IVD	-	13q12.11	●	20 (200 µl)
Z-2095-50/-200	ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe C€ IVD	D18Z1	13q12.11/18p11.1/21q22.13-q22.2	●/●/●	5/20 (50/200 µl)
Z-2164-200	ZytoLight SPEC 13/21 Dual Color Probe C€ IVD	-	13q12.11/21q22.13-q22.2	●/●	20 (200 µl)
Z-2006-200	ZytoLight CEN 17 Probe C€ IVD	D17Z1	17p11.1-q11.1	●	20 (200 µl)
Z-2007-200	ZytoLight CEN 18 Probe C€ IVD	D18Z1	18p11.1-q11.1	●	20 (200 µl)
Z-2163-200	ZytoLight SPEC 18/CEN X/Y Triple Color Probe C€ IVD	DXZ1/DYZ3	18q21.31-q21.32/Xp11.1-q11.1/Yp11.1-q11.1	●/●/●	20 (200 µl)
Z-2086-200	ZytoLight SPEC 21q22 Probe C€ IVD	-	21q22.13-q22.2	●	20 (200 µl)
Z-2180-200	ZytoLight SPEC 21/CEN X/Yq12 Triple Color Probe C€ IVD	DXZ1/ III DYZ1	21q22.13-q22.2/Xp11.1-q11.1/Yq12	●/●/●	20 (200 µl)
Z-2008-200	ZytoLight CEN X Probe C€ IVD	DXZ1	Xp11.1-q11.1	●	20 (200 µl)
Z-2010-200	ZytoLight CEN Yq12 Probe C€ IVD	III DYZ1	Yq12	●	20 (200 µl)
Z-2123-200	ZytoLight CEN Y (DYZ3) Probe C€ IVD	DYZ3	Yp11.1-q11.1	●	20 (200 µl)
Z-2016-50/-200	ZytoLight CEN X/Yq12 Dual Color Probe C€ IVD	DXZ1/III DYZ1	Xp11.1-q11.1/Yq12	●/●	5/20 (50/200 µl)
Z-2120-200	ZytoLight CEN X/Y Dual Color Probe C€ IVD	DXZ1/ DYZ3	Xp11.1-q11.1/Yp11.1-q11.1	●/●	20 (200 µl)
Related Products					
Z-2279-20	ZytoLight Aneuploidy Panel 18/X/Y and 13/21 C€ IVD Incl. ZytoLight SPEC 18/CEN X/Y Triple Color Probe, 0.2 ml (Z-2163-200); ZytoLight SPEC 13/21 Dual Color Probe, 0.2 ml (Z-2164-200)		✓		20
Z-2104-5	ZytoLight Aneuploidy Panel X/Y and 13/18/21 C€ IVD Incl. ZytoLight CEN X/Yq12 Dual Color Probe, 0.05 ml (Z-2016-50); ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe, 0.05 ml (Z-2095-50)		✓		5
Z-2104-20	ZytoLight Aneuploidy Panel X/Y and 13/18/21 C€ IVD Incl. ZytoLight CEN X/Yq12 Dual Color Probe, 0.2 ml (Z-2016-200); ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe, 0.2 ml (Z-2095-200)		✓		20
Z-2028-5	ZytoLight FISH-Tissue Implementation Kit C€ IVD Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 150 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml				5
Z-2028-20	ZytoLight FISH-Tissue Implementation Kit C€ IVD Incl. Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 500 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml				20
Z-2099-20	ZytoLight FISH-Cytology Implementation Kit C€ IVD Incl. Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x MgCl ₂ , 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/DuraTect-Solution, 0.8 ml				20



Pragostem s.r.o.

Provozovna: 17.listopadu 2753, 530 02 Pardubice

E-mail: info@pragostem.cz

E-mail: objednavky@pragostem.cz

www.pragostem.cz